

What's Inside?



How Does It Work?

info@PediEd.com

888-280-PEDS(7337)

PediEd.com



30 life-size foldable, durable, realistic "Mani-Kids" These pediatric patients range from 6-months to 12-years-old. Injuries including everything from minor scrapes and bruises to severe burns and amputations.



5 "Mega Mani-Kids" AKA adults



3 "Mani-Mutts"

These are a great opportunity to talk with your students about specific protocols regarding animals.



4 Triage Tarps



30 Pedi-Wheels and 12 clipboards



EF4 Tornado at the Drive-In Airplane Crashes on Takeoff Elementary/Middle School Shooting 30-Car Pileup On the Highway Fire at a Daycare/Afterschool Program Catastrophic Rollercoaster Collapse

School Bus Grashes Off Bridge Into River Below Truck Plows Through Pizza Palace Playland Sports Complex Bombing Flying Bouncy Houses at Easter Picnic Massive 7.4 Earthquake Devastates Small Town Gas Leak with Explosion at Build-A-Bunny

12 unique, real life, disaster scenarios!

EF4 Tornado at the Drive-In



Inside <u>each</u> of the 12 clear color-coded disaster scenario folders you will find... Scenario specific triage cards, removable moulage, answer keys, and sound effects!



Scenario



It's the much-anticipated opening weekend at the drive-in theater just outside of form. During the first movie, The Wizard of O.z. community-based emergency alorm sires started going off as an explosion from a large propane tank and fire in the concrete concession stand where many families were taking retuge. Others, still in their cas, were filed and thrown around the area. Some of the smaller vehicles were thrown into trees and ditches. There are several injured adults: however, the majority of the vicinities are young children. *1 Hour Of Scenario Specific Sound Effects!*





The front of the educator booklet. You will find the scenario and the QR code for 1-hour of scenario specific sound effects!



Inside the educator booklet, you will find which Mani-Kids are used for the scenario, what temporary moulage pieces they need, and what their triage color is for that specific scenario.

EF4 Tornado at the Drive-In

Student Answer Form—Tornado MCI



ige Color? ige Color? ige Color? ige Color? ige Color? ige Color? ige Color? ige Color?	Question 2: Question 2: Question 2: Question 2: Question 2: Question 2:	Question 3: Question 3: Question 3: Question 3: Question 3:		Mani-Kid #2 Mani-Kid #3 Mani-Kid #4 Mani-Kid #5	Triage Color? Triage Color? Triage Color? Triage Color?	Red Green	Question 2: Question 2: Question 2:	Red F: 49 mcg -OR- K: 1.6-9 mg -OR- M: 0.4-0.9 mg* Purple	Question 3: Question 3: Question 3:	BP 66-126 3.5-4.0 uncuffed ETT -OR- 3.0-3.5 cuffed ETT 10-11 kg
ige Color? ige Color? ige Color? ige Color? ige Color? ige Color?	Question 2: Question 2: Question 2:	Question 3: Question 3:		Mani-Kid #4	Triage Color?	Green	Question 2:	Purple	Question 3:	-OR- 3.0-3.5 cuffed ETT
ige Color? ige Color? ige Color? ige Color?	Question 2: Question 2:	Question 3:			and the second sec			Contract of the second		10-11 kg
ige Color? ige Color? ige Color?	Question 2:			Mani-Kid #5	Triage Color?	NAME AND ADDRESS OF				
ige Color? ige Color?		Question 3:			muge color:	Yellow	Question 2:	RR 20-30	Question 3:	-OR- M-0.6-1.4 mg*
ige Color?	Question 2:		\leftarrow	Mani-Kid #6	Triage Color?	Red	Question 2:	240-290 mL Normal Saline or LR	Question 3:	-OR- M: 0.6-1.4 mg*
	and a structure and	Question 3:	• P	Mani-Kid #7	Triage Color?	Green	Question 2:	HR 65-140	Question 3:	F:8-18 mcg -OR- K:3-18 m -OR- M: 0.8-1.8 mg
	Question 2:	Question 3:		Mani-Kid #8	Triage Color?	Red	Question 2:	5.0 uncuffed ETT -OR- 4.5 cuffed ETT	Question 3:	8-10 FR
ige Color?	Question 2:	Question 3:		Mani-Kid #9	Triage Color?	Yellow	Question 2:	White	Question 3:	F: 8-18 mcg -OR- K: 3-18 m -OR- M: 0.8-1.8 mg
ge Color?	Question 2:	Question 3:		Mani-Kid #10	Triage Color?	Red	Question 2:	Absolutely, YES	Question 3:	*NASAL F: 15-36 mcg -OR- K: 15-90 r
ge Color?	Question 2:	Question 3:	\leftarrow	Mani-Kid #11	Triage Color?	Black	Question 2:	King Sze 2.5 -OR- LMA Sze 2-2.5 -OR- i-gel Sze 2	Question 3:	1.9-2.3 mL EPINEPHrine
ge Color?	Question 2:	Question 3:		Mani-Kid #12	Triage Color?	Yellow	Question 2:	Blue	Question 3:	19-23 kg
ge Color?	Question 2:	Question 3:		Mani-Kid #13	Triage Color?	Red	Question 2:	5.5 uncuffed ETT -OR- 5.0 cuffed ETT	Question 3:	15-16.5 cm
ge Color?	Question 2:	Question 3:		Mani-Kid #14	Triage Color?	Red	Question 2:	380-460 mL Normal Saline or LR	Question 3:	P: 10-23 mcg -OR- K: 3.8-23 -OR- M: 1-2.3 mg*
ge Color?	Question 2:	Question 3:		Mani-Kid #15	Triage Color?	Red	Question 2:		Question 3:	11.4-27.6 mg Rocuronium
ge Color?	Question 2:	Question 3:	\leftarrow	Mani-Kid #16	Triage Color?	Black	Question 2:	19-23 kg	Question 3:	1.9-2.3 mL EPINEPHripp
ge Color?	Question 2:	Question 3:	• • • •	Mani-Kid #17	Triage Color?	Green	Question 2:	HR 60-130	Question 3:	RR 12-22
ge Color?	Question 2:	Question 3:		Mani-Kid #18	Triage Color?	Red	Question 2:	480-580 mL Normal Saline or LR	Question 3:	F: 12-29 mcg -OR- K: 48-29 i -OR- M: 1.2-2.9 mg
ge Color?	Question 2:	Question 3:		Mani-Kid #19	Triage Color?	Yellow	Question 2:	* NASAL V: 24-58 mg -OR- A: 24-29 mg	Question 3:	480-580 mL Normal Saline or LR
ge Color?	Question 2:	Question 3:		Mani-Kid #20	and the second se		Question 2:	24-29 kg	Question 3:	-OR- M: 1.2-29 mg*
ge Color?	Question 2:	Question 3:	•	Mani-Kid #21	-		Question 2:	HR 60-110	Question 3:	BP 88-120
ge Color?	Question 2:	Question 3:		Mani-Kid #22	Triage Color?	Red	Question 2:	F: 15-36 mcg -OR- K: 6-36 mg	Question 3:	18-43.2 mg Rocuronium
ge Color?	Question 2:	Question 3:		Mani-Kid #23	Triage Color?	Red	Question 2:	* MATAI	Question 3:	600-720 mL Normal Saline or LR
ge Color?	Question 2:	Question 3:		Mani-Kid #24	a second and a second second	Yellow	Question 2:		Question 3:	30-36 kg
ge Color?	Question 2:	Question 3:		Mani-Kid #25	and the second sec		Question 2:	King Size 2.5-3 -OR- LMA Size 3	Question 3:	3-3.6 mL
ge Color?	Question 2:	Question 3:		Mani-Kid #26	and the second se		Question 2:	BP 88-120	Question 3:	*NASAL F: 30-80 mcg -OR- K: 30-200 r
ge Color?	Question 2:	Question 3:		Mani-Kid #27	Triage Color?	Green	Question 2:	HR 60-110	Question 3:	RR 12-22
ge Color?	Question 2:	Question 3:	\leftarrow	Mani-Kid #28			Question 2:		Question 3:	F: 20-50 mcg -OR- K: 8-50 m -OR- M: 2-5 mg*
ge Color?	Question 2:	Question 3:		Mani-Kid #29			Question 2:	6.5-7.0 cuffed ETT	Question 3:	20-22 cm
ige Color?	Question 2:	Question 3:		Mani-Kid #30			Question 2:		Question 3:	12-14 FR
ge Color?	Question 2:	Question 3:		Mega Mani-Kid #31	A REAL PROPERTY AND A REAL PROPERTY.		Question 2:		Question 3:	-25-100 mcg -OR- K: 20-100 -OR- M: 2-10 mg (ADULT)
ge Color?	Question 2:	Question 3:		Mega Mani-Kid #32	Triage Color?	Yellow	Question 2:	F: 25-100 mcg -OR- K: 20-100 mg		500-1000 mL Normal Saline or (ADULT)
ge Color?	Question 2:	Question 3:		and the second s	and the second sec		Question 2:	HR 60-100	Question 3:	RR 12-20
ge Color?	Question 2:	Question 3:		International Content of the International	And the second second		Question 2:	* (ADULT) NASAL F: 50-100 mcg -OR- 1: 25-100 mg		-25-100 mcg -OR- K: 20-100 -OR- M: 2-10 mg (ADULT)
	ge Color? ge Color?	ge Color? Question 2: ge Color? Qu	ge Color? Question 2: Question 3: ge Color? Question 2: Que	ge Color? Question 2: Question 3: ge Color? Question 2: Que	ge Color? Question 2: Question 3: Mani-Kid #12 ge Color? Question 2: Question 3: Mani-Kid #13 ge Color? Question 2: Question 3: Mani-Kid #14 ge Color? Question 2: Question 3: Mani-Kid #15 ge Color? Question 2: Question 3: Mani-Kid #16 ge Color? Question 2: Question 3: Mani-Kid #17 ge Color? Question 2: Question 3: Mani-Kid #18 Mani-Kid #17 Mani-Kid #18 Mani-Kid #18 Mani-Kid #20 Mani-Kid #20 Mani-Kid #21 Mani-Kid #20 Question 3: Mani-Kid #23 ge Color? Question 2: Question 3: Mani-Kid #24 ge Color? Question 2: Question 3: Mani-Kid #24 Mani-Kid #24 Mani-Kid #24 Mani-Kid #25 Mani-Kid #26 Mani-Kid #24 Question 3: Mani-Kid	ge Color? Question 2: Question 3: Mani-Kid #12 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #13 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #15 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #16 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #12 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #20 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #24 Triage Color? ge Color? Question 2: Question 3: Mani-Kid #24 Triage Color?	ge Color? Question 2: Question 3: ge Color? Question 2: Que	ge Color? Question 2: Question 3: Mani-Kid #12 Triage Color? Yellow Question 2: ge Color? Question 2: Question 3: Mani-Kid #13 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #16 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #16 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? Red Question 2: ge Color? Question 2: Question 3: Mani-Kid #21 Triage Color? Red Question 2:<	ge Color? Question 2: Question 3: Mani-Kid #12 Triage Color? Yellow Question 2: Blue ge Color? Question 2: Question 3: Mani-Kid #13 Triage Color? Red Question 2: Blue ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? Red Question 2: Horn Science vic Red ge Color? Question 2: Question 3: Mani-Kid #14 Triage Color? Red Question 2: Horn Science vic Red ge Color? Question 2: Question 3: Mani-Kid #15 Triage Color? Red Question 2: HP-23 kg Mani-Kid #15 Triage Color? Red Question 2: HR 60-130 ge Color? Question 2: Question 3: Mani-Kid #17 Triage Color? Red Question 2: HR 60-130 ge Color? Question 2: Question 3: Mani-Kid #19 Triage Color? Red Question 2: HR 60-110 ge Color? Question 2: Question 3: Mani-Kid #23 Triage Color? Red Question 2: HR 60-110 ge Color? Question 2:	ge Color? Question 2: Question 3: Mani-Kid #12 Triage Color? Yellow Question 2: Blue Question 3: ge Color? Question 2: Question 3: Mani-Kid #13 Triage Color? Red Question 2: Question 3: Question 3:

On the last page of all the educator booklets, you will find... The blank <u>scenario specific</u> "Student" answer sheet to make photocopies, and the "Educator" answer key for easy grading.



Inside the scenario folder, you will find a pack of triage cards. Each Mani-Kid has their own card for each MCI. Students will reference the information on the cards, as well as the visual moulage on the patient to determine the correct triage color.



Triage card breakdown

Note: The layout is the same for every card, in every one of the 12 disaster scenarios.



The Mani-Kid number that is in the top right corner of the board also corresponds to the triage card number.



After your students have triaged all the patients, there are two additional pediatric questions on the back of the card.

> *Note: These <u>are not</u> secondary triage questions... They are intended for additional pediatric practice only.*



The two additional questions change depending on the scenario and the injuries of the patient.

Note: The questions <u>do not</u> indicate the next steps of care, i.e., if the kid is triaged to green and question #3 is "what size ET tube for this patient" that does not mean you should be putting an ET tube in a healthy kid. Again, these questions are there to give students the opportunity to use their cheat sheets and ultimately practice <u>finding</u> the information.



Depending on the scenario, certain Mani-Kids may require temporary moulage. The scenario specific educator manual will show which kids need what moulage and where.



<u>Card Information</u> Child is walking around with blood on her right hand/wrist

Temporary moulage pieces are clear, sticky, plastic overlays to be placed on the Mani-Kids. The moulage matches the triage card patient information. Note: Inside each scenario folder, there is <u>all</u> the moulage for <u>all</u> the Mani-Kids for <u>that specific scenario</u>.



Running one disaster scenario is as easy as 1-2-3-4

- 1- Pick one MCI scenario (example, The Tornado MCI)
- 2- Match the triage cards with their Mani-Kids
- 3- Add additional temporary moulage if needed
- 4- Put the patients around the room, ambulance bay, parking lot, etc. That's It! Your ready to train!



Each student will need: A clipboard, scenario specific answer sheet, (i.e., The Tornado MCI), Pedi-Wheel and a pen! That's it!

Note: Students are also encouraged to use apps on their phone or other cheat sheets as needed to determine correct triage color.



They are now ready to practice their triaging skills and gain additional pediatric knowledge! They also get to practice moving patients to appropriate tarps or areas using team work, and ultimately through training on a more regular basis, be more prepared for a real-life disaster.



But it doesn't end there! Educators, take things up a notch by turning off the lights, adding a fog machine, flashing lights, the included sound effects and making your students use a flashlight! The more intense the training... The stronger they will become!



Verify all patients were triaged correctly using the triage answer key inside the <u>scenario specific</u> educator manual.



Credit hours? Have your students log in to the online component of Disasters In Seconds, fill in their answers at home and get 1 credit hour for each disaster scenario!





With 12 unique disaster scenarios, you could run a completely hands-on MCI training day once a month for an entire year!



There is a time and place for real life actors, makeup artists and gooey, drippy moulage when it comes to disaster training... because it's cool!! But trying to coordinate with schedules, paying for overtime, and getting patients to act out their actual injuries can be a LOT of work. Disasters In Seconds alleviates all of these issues, takes one person about 20-30 minutes to set up, and at the end of the day all fits back in one bag. DISASTER TRAINING DONE!!



Disaster Drills & Triaging Tiny Traumas info@PediEd.com 888-280-PEDS(7337) PediEd.com

✓ 35 Patients

I2 MCI Scenarios

I Disaster Training System
You'll Never Forget!



